



ANAEROBIC JARS LAJ-B1 SERIES

ANAEROBIC JARS LAJ-B1 SERIES

ANAEROBIC JAR LAJ-B10

Anaerobic Jar LAJ-B10 is a pumping type anaerobic jar with 7 L of capacity and 3/25 R petri dish holder. Featuring PMMA transparent material with extremely smooth surface along with excellent airtight performance and avoidance of any undue leakage. Designed with conventional O-rings and screw clamp to ensure a more positive seal. Equipped with snap-shut gas supply coupling for easy connecting.

ANAEROBIC JAR LAJ-B11

Anaerobic Jar LAJ-B11 is a hermetic type anaerobic jar with 3/15 R petri dish holder. Featuring PMMA transparent material with extremely smooth surface along with excellent airtight performance and avoidance of any undue leakage. Designed with conventional O-rings and screw clamp to ensure a more positive seal. Equipped with snap-shut gas supply coupling for easy connecting.

FEATURES

- Transparent lid and jar with O-ring inbetween
- Easy to clean
- Lightweight design with convenient portability
- Optional racks providing holders for comfortable operation of the anaerobe systems
- Durable stainless steel clamp with power electrostatic painting
- Another side tube for evacuation and introduction of gases
- Shock and impact resistance
- Durable and convenient for experimental management

APPLICATIONS

Anaerobic Jars are used for anaerobic incubation (anaerobiosis) of bacterial cultures for research and testing across the fields of microbiology, biotechnology, pathology, virology and CDC.

SPECIFICATION

| Model | LAJ-B10 | LAJ-B11 |
|--|--|--|
| Type | Pumping type | Hermetic type |
| Time taken (to attain anaerobic condition) | 3 to 5 mins | 4 hours |
| O-ring | Yes | Yes |
| Contents | 3 stacks of 08 dishes (ϕ 9-10cm) | 3 stacks of 08 dishes (ϕ 9-10cm) |
| Petridish holder | 3/25R | 3/15R |
| Material | PMMA material | PMMA material |
| Clamp | Colored stainless steel clamp | Colored stainless steel clamp |
| Quick snap-shut coupling | Yes | Yes |
| Packaging dimensions (L×W×H) | 370×370×370 mm | 370×370×370 mm |
| Gross weight | 6.3 kg | 6 kg |