

# Manual Dual Door Flammable Cabinets

## LFSC-A1 Series



### 114 L Flammable Storage Cabinet LFSC-A10

---

114 L Flammable Storage Cabinet LFSC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- ◆ Capacity - 114 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ One adjustable shelf
- ◆ Adjustable spill retaining shelves

### 170 L Flammable Storage Cabinet LFSC-A11

---

170 L Flammable Storage Cabinet LFSC-A11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- ◆ Capacity - 170 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ Two adjustable shelves
- ◆ Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specification

Model No	LFSC-A10	LFSC-A11
Capacity	114 L	170 L
Door type	Dual door manual	Dual door manual
Adjustable shelf	1	2
Shelves loading capacity	50 kgs	100 kgs
Dimension ( H x W x D )	1120 x 1090 x 460 mm	1650 x 1090 x 460 mm
Weight	107 kgs	148 kgs

## 204 L Flammable Storage Cabinet LFSC-A12

204 L Flammable Storage Cabinet LFSC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

### Features

- ◆ Capacity - 204 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ One adjustable shelf
- ◆ Adjustable spill retaining shelves

## 227 L Flammable Storage Cabinet LFSC-A13

227 L Flammable Storage Cabinet LFSC-A13 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

### Features

- ◆ Capacity - 227 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ Two adjustable shelves
- ◆ Adjustable spill retaining shelves

### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

### Specification

Model No	LFSC-A12	LFSC-A13
Capacity	204 L	227 L
Door type	Dual door manual	Dual door manual
Adjustable shelf	1	2
Shelves loading capacity	100 kgs	100 kgs
Dimension ( H x W x D )	1650 x 860 x 860 mm	1650 x 860 x 860 mm
Weight	202 kgs	173 kgs

### 340 L Flammable Storage Cabinet LFSC-A14

---

340 L Flammable Storage Cabinet LFSC-A14 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- ◆ Capacity - 340 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ Two adjustable shelves
- ◆ Adjustable spill retaining shelves

### 415 L Flammable Storage Cabinet LFSC-A15

---

415 L Flammable Storage Cabinet LFSC-A15 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- ◆ Capacity - 415 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ Two adjustable shelves
- ◆ Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specification

Model No	LFSC-A14	LFSC-A15
Capacity	340 L	415 L
Door type	Dual door manual	Dual door manual
Adjustable shelf	2	2
Shelves loading capacity	100 kgs	100 kgs
Dimension ( H x W x D )	1650 x 1090 x 860 mm	1650 x 1500 x 860 mm
Weight	186 kgs	293 kgs

### 434 L Flammable Storage Cabinet LFSC-A16

434 L Flammable Storage Cabinet LFSC-A16 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity

#### Features

- ◆ Capacity - 434 L
- ◆ Dual walled manual
- ◆ 38 mm insulating air space
- ◆ 180° fully opening door
- ◆ Three point latch and manual lock
- ◆ Three adjustable shelves
- ◆ Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Model No	LFSC-A16
Capacity	434 L
Door type	Dual door manual
Adjustable shelves	3
Shelves loading capacity	100 kgs
Dimension ( H x W x D )	1650 x 1500 x 860 mm
Weight	290 kgs