

Manual Dual Door
Weak Acid and
Alkali Cabinet
Storage Cabinet
LWAC-A1 Series



Manual Dual Door Weak Acid and Alkali Chemical Storage Cabinet

114 L Weak Acid and Alkali Cabinet LWAC-A10

114 L Weak Acid and Alkali Cabinet LWAC-A10 is a dual walled Weak Acid and Alkali resistance cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 114 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-A10
Capacity	114 L
Door type	Dual door manual
Adjustable shelf	1
Dimension (W x D x H)	1090 x 460 x 1120 mm
Weight	117 kgs

Manual Dual Door Weak Acid and Alkali Chemical Storage Cabinet

170 L Weak Acid and Alkali Cabinet LWAC-A11

170 L Weak Acid and Alkali Cabinet LWAC-A11 is a dual walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 170 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-A11
Capacity	170 L
Door type	Dual door manual
Adjustable shelf	2
Dimension (W x D x H)	1090 x 460 x 1650 mm
Weight	155 kgs

Manual Dual Door Weak Acid and Alkali Chemical Storage Cabinet

227 L Weak Acid and Alkali Cabinet LWAC-A12

227 L Weak Acid and Alkali Cabinet LWAC-A12 is a dual walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 227 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-A12
Capacity	227 L
Door type	Dual door manual
Adjustable shelf	2
Dimension (W x D x H)	1650 x 860 x 860 mm
Weight	173 kgs

Manual Dual Door Weak Acid and Alkali Chemical Storage Cabinet

340 L Weak Acid and Alkali Cabinet LWAC-A13

340 L Weak Acid and Alkali Cabinet LWAC-A13 is a dual walled manual Weak Acid and Alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 340 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications

Model No	LWAC-A13
Capacity	340 L
Door type	Dual door manual
Adjustable shelves	2
Dimension (H x W x D)	1650 x 1090 x 860 mm
Weight	186 kgs