

Manual Single Door

Flammable Storage Cabinet

LFSC-B Series



15 L Flammable Storage Cabinet LFSC-B10

15 L Flammable Storage Cabinet LFSC-B10 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 15 L
 - Dual walled manual
 - 38 mm insulating air space
 - 180° fully opening door
 - Three point latch and manual lock
 - One adjustable shelf
 - Adjustable spill retaining shelves
-

38 L Flammable Storage Cabinet LFSC-B11

38 L Flammable Storage Cabinet LFSC-B11 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 38 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specification

Model No	LFSC-B10	LFSC-B11
Capacity	15 L	38 L
Door type	Single door manual	Single door manual
Adjustable shelf	1	1
Shelves loading capacity	50 kgs	50 kgs
Dimension (H x W x D)	560 x 430 x 430 mm	640 x 590 x 600 mm
Weight	30.5 kgs	56 kgs

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

45 L Flammable Storage Cabinet LFSC-B12

45 L Flammable Storage Cabinet LFSC-B12 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 45 L
- Dual walled manual
- 38 mm insulating air space
- 1800 fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

83 L Flammable Storage Cabinet LFSC-B13

83 L Flammable Storage Cabinet LFSC-B13 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 83 L
- Dual walled manual
- 38 mm insulating air space
- 1800 fully opening door
- Three point latch and manual lock
- Three adjustable shelves
- Adjustable spill retaining shelves

Specification

Model No	LFSC-B12	LFSC-B13
Capacity	45 L	83 L
Door type	Single door manual	Single door manual
Adjustable shelf	1	3
Shelves loading capacity	50 kgs	50 kgs
Dimension (H x W x D)	890 x 590 x 460 mm	1650 x 600 x 460 mm
Weight	59 kgs	99 kgs

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

204 L Flammable Storage Cabinet LFSC-B14

204 L Flammable Storage Cabinet LFSC-B14 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 204 L
 - Dual walled manual
 - 38 mm insulating air space
 - 180° fully opening door
 - Three point latch and manual lock
 - Three adjustable shelves
 - Adjustable spill retaining shelves
-

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specification

Model No	LFSC-B14
Capacity	204 L
Door type	Single door manual
Adjustable shelves	3
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 600 x 870 mm
Weight	150 kgs