

# PHARMACY REFRIGERATOR

## LPRF-B1 SERIES



### **Pharmacy Refrigerator LPRF-B10**

Pharmacy Refrigerator LPRF-B10 is a digital microcomputer-controlled refrigerator with holding capacity of 660 L. Features branded compressor and fan with 8 to 20°C of temperature range. Ergonomic design with in-built lamp and safety lock, has noiseless and user-friendly operation mode. With auto audio-visual signal alarm protection and real-time temperature display, this unit offers stable and reliable performance for storing different biological products.

### **Pharmacy Refrigerator LPRF-B11**

Pharmacy Refrigerator LPRF-B11 is a digital microcomputer-controlled refrigerator with holding capacity of 868 L. Features branded compressor and fan with 8 to 20°C of temperature range. Ergonomic design with in-built lamp and safety lock, has noiseless and user-friendly operation mode. With auto audio-visual signal alarm protection and real-time temperature display, this unit offers stable and reliable performance for storing different biological products.

### **Pharmacy Refrigerator LPRF-B12**

Pharmacy Refrigerator LPRF-B12 is a digital microcomputer-controlled refrigerator with holding capacity of 960 L. Features branded compressor and fan with 8 to 20°C of temperature range. Ergonomic design with in-built lamp and safety lock, has noiseless and user-friendly operation mode. With auto audio-visual signal alarm protection and real-time temperature display, this unit offers stable and reliable performance for storing different biological products.

### **Pharmacy Refrigerator LPRF-B13**

Pharmacy Refrigerator LPRF-B13 is a digital microcomputer-controlled refrigerator with holding capacity of 1008 L. Features branded compressor and fan with 8 to 20°C of temperature range. Ergonomic design with in-built lamp and safety lock, has noiseless and user-friendly operation mode. With auto audio-visual signal alarm protection and real-time temperature display, this unit offers stable and reliable performance for storing different biological products.

### **Pharmacy Refrigerator LPRF-B14**

Pharmacy Refrigerator LPRF-B14 is a digital microcomputer-controlled refrigerator with holding capacity of 1360 L. Features branded compressor and fan with 8 to 20°C of temperature range.

Ergonomic design with in-built lamp and safety lock, has noiseless and user-friendly operation mode. With auto audio-visual signal alarm protection and real-time temperature display, this unit offers stable and reliable performance for storing different biological products.

## Features

- ❑ A digital microcomputer-controlled refrigerator with large capacity for storage
- ❑ Branded compressor and fan with 8 to 20°C of temperature range, effective refrigeration
- ❑ Ergonomic design with in-built lamp and safety lock, convenient and safe storage
- ❑ Noiseless and user-friendly operation mode, efficient and reliable operation
- ❑ Auto audio-visual signal alarm protection, ensures safe and secure operation
- ❑ Fluorine-free refrigerant, environment-friendly and efficient refrigeration
- ❑ Stable and reliable performance for storing different biological products

## Application

Pharmaceutical Refrigerator is used to store various drugs, reagents, vaccines and biological products across medical institutions, pharmaceutical factories, hospitals, health clinics, scientific research institutes, epidemic prevention centers, and other industrial institutions etc.

## Specifications (LPRF-B10, LPRF-B11 and LPRF-B12)

Model	LPRF-B10	LPRF-B11	LPRF-B12
Cabinet	Upright		
Capacity	660 L	868 L	960 L
Temperature range	8 to 20°C		
Compressor	1 pc		
Refrigerant	Fluorine-free		
Defrost mode	Automatic		
Alarm system	High/low temperature, Power off, Door ajar, Sensor fault etc.		
Controller	Microcomputer digital controller		
Power supply	AC 187 to 242 V		
Power	380 W	450 W	

## PHARMACY REFRIGERATOR LPRF-B1 SERIES

Internal dimension (W×D×H)	1112×560×1232 mm	1100×634×1325 mm	1080×610×1380 mm
External dimension (W×D×H)	1200×650×1895 mm	1220×755×1885 mm	1200×745×1870 mm
Net weight	129 kg	178 kg	141 kg

### Specifications (LPRF-B13 and LPRF-B14)

Model	LPRF-B13	LPRF-B14
Cabinet	Upright	
Capacity	1008 L	1360 L
Temperature range	8 to 20°C	
Compressor	1 pc	
Refrigerant	Fluorine-free	
Defrost mode	Automatic	
Alarm system	High/low temperature, Power off, Door ajar, Sensor fault etc.	
Controller	Microcomputer digital controller	
Power supply	AC 187 to 242 V	
Power	480 W	650 W
Internal dimension (W×D×H)	1100×734×1325 mm	1680×610×1380 mm
External dimension (W×D×H)	1220×855×1885 mm	1800×745×1870 mm
Net weight	194 kg	158 kg